



Synergent Biochem Inc.

LIQUID LIPID CONTROL

For *in vitro* diagnostic use

Cat. No. 915, Levels 1 & 2

INTENDED USE

The Synergent Biochem, Inc. bi-level Lipid Control is a ready-to-use, assayed chemistry control material for use in monitoring the accuracy and precision of the quantitative lipid determinations.

SUMMARY AND PRINCIPLE

Synergent Liquid Lipid control is prepared from fresh human sera enriched with human lipid fractions. This assayed control includes assigned values for total cholesterol, direct HDL and LDL cholesterol and triglycerides. The use of assayed quality control materials in clinical chemistry procedures is necessary in maintaining good laboratory practices.

PRECAUTIONS

Human source material: Treat as potentially infectious. Each donor unit used in the preparation of this product has been tested by FDA accepted methods and found non-reactive for hepatitis B surface antigen (HBsAg), for antibody against human immune-deficiency virus (anti-HIV-1/2), and for antibody against hepatitis C virus (anti HCV). In spite of negative results at donor level, this control should be handled as potentially infectious material in accordance with good laboratory practices (GLP).

PROCEDURE

1. Remove control from refrigerator and allow warming to room temperature for 15 to 30 minutes.
2. Gently swirl the contents until homogenous. Do not shake vial to avoid the formation of foam.
3. Recap vial tightly after use. Store opened vial in refrigerator at 2° to 8° C

MATERIALS PROVIDED

Liquid Lipid Control - 6 vials, 3mL/vial.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Reagents and supplemental material required for the measurement of listed constituents.
2. General laboratory equipment
3. Suitable analyzer system.

PRODUCT STABILITY

The Synergent Liquid Lipid Control is stable until the expiration date when stored at +2 °C to +8 °C. After opening the vial, the control is stable for 30 days when stored at +2 °C to +8 °C.

LIMITATIONS OF PROCEDURE

Erroneous results may occur from:

1. Contamination of control. Discard vial if particulates or other signs of microbial contamination are observed.
2. Technique errors associated with reagent preparation and the assay procedure, using this product as a calibrator or using past the expiration date.

VALUE ASSIGNMENT

The mean values were derived from replicate analyses of vials representative of the entire lot. The expected range of the mean is provided as a reference for the laboratory until it has established its own mean and standard deviation. The indicated mean and its assigned range should serve as guidelines in assessing the performance of each test method. Measurements using other reagents and instrument systems may give different results.

ASSAY VALUES

| Analyte / Instrument | Unit | Level 1 | | Level 2 | |
|---------------------------|-------|------------------------------------|---------|------------------------------------|---------|
| | | Lot: 1105702E Expiry: 02 / 2013 | | Lot: 1105703E Expiry: 02 / 2013 | |
| | | Mean | Range | Mean | Range |
| Cholesterol, Total | | | | | |
| ChemWell | mg/dL | 193 | 154-232 | 300 | 240-360 |
| Roche Modular | mg/dL | 215 | 172-258 | 303 | 242-364 |
| Triglycerides | | | | | |
| ChemWell | mg/dL | 145 | 116-174 | 394 | 315-473 |
| Roche Modular | mg/dL | 134 | 107-161 | 372 | 298-446 |
| Cholesterol, HDL | | | | | |
| ChemWell | mg/dL | 66 | 53-79 | 41 | 33-49 |
| Roche Modular | mg/dL | 73 | 58-88 | 40 | 32-48 |
| Cholesterol, LDL | | | | | |
| ChemWell | mg/dL | 107 | 86-128 | 191 | 153-229 |
| Roche Modular | mg/dL | 122 | 98-146 | 204 | 163-245 |
| Apolipoprotein A1 | | | | | |
| Beckman Immage | mg/dL | \$ | | \$ | |
| Roche Modular | mg/dL | 143 | 114-172 | 90 | 72-108 |
| Apolipoprotein B | | | | | |
| Beckman Immage | mg/dL | \$ | | \$ | |
| Roche Modular | mg/dL | 98 | 78-118 | 168 | 134-202 |

§ The data required to establish the means and acceptable ranges for this assay were not obtained due to limited assignment participation.

TECHNICAL ASSISTANCE

For technical assistance and ordering information contact Synergent Biochem Inc. at 562-809-3389.

CATALOG NUMBERS:

915 6 X 3mL BI-LEVEL
915-1 6 X 3mL LEVEL 1
915-2 6 X 3mL LEVEL 2

Rev10/07